

Pharmacy Refrigerator and Freezer

HP series fridges designs for demands for refrigerator chamber and freezer chamber separately in one cabinet. It can be used for storage pharmacy, vaccine, medicine, blood products, biology products, and so on for hospital, blood-bank, center for Disease Prevention, clinic and scientific research institution, etc.





1,Prices Temperature control

Refrigerator Chamber: 2° C ~ 8° C;

Freezer Chamber:- 10° C ~ -26 $^{\circ}$ C & - 10° C ~ -40 $^{\circ}$ C;

Controlled by micro-computer; temperature for refrigerator and freezer can be digital display separately;

Forced air cooling system, the refrigeration function can be stop separately.

2,Safety Control Systems

Light and sound alarms:high or low temperature,temperature sensor failure,door open(For refrigerator);

Protecting functions such as power-on delay, off-cycle protection.



3,Cooling system

Select famous brand compressor and introduced Fan Motor,rational optimized evaporator and the new wind structure, ensure uniform temperature stability in refrigerator;

Fluoride-free foaming, fluoride-free refrigerant, environment-friendly.

4, Humanized Design

LCD display(for PF210 and PF265),LED display(for PF-288 and PF-288A);

With shelves(PF210/PF265 is glass shelves) for refrigerator chamber, with transparent drawer for freezer chamber;

Safety door lock, preventing the door opened by chance;

Universal caster with locking function design at the bottom, free to moved and easy to fix; Wide range of voltage design, range from 198~242 volts, 50HZ/60HZ(Optional).

Main Technical Parameters:

Model	Capacit y	Temn	External Size	Inside Size(W*D*H/MM)			Shelves
				Refrigerator	Free zer	Climate	/Drawe r
HP-PF210	210L	2 ~ 8 °C; -10 ~ -26 °C	555*580*1670	428*430*830	380*420*560	N	3/4
HP-PF270	2701	2 ~ 8 °C; -10 ~ -26 °C	600*590*1850	499*520*690	460*469*650	N	3/4
HP-PF290	70M	2 ~ 8 °C; -10 ~ -40 °C	740*780*1850	590*550*630	510*490*420	N	3/3
HP-PF290A	1290L	2 ~ 8 °C; -10 ~ -40 °C	740*780*1850	590*550*630	510*490*420	N	3/3

Note:More information for our products and service, please turn to our official website or contact us directly.