

DW7070H Laboratory Opener



MODEL NO.

DW7010H Laboratory Opener

Application:

DW7070H Laboratory Opener is used to make opening and scotching raw cotton to be monofilament and remove the impurities. It's the ideal system for teaching, researching, etc.

Features:

- High-efficiency imported frequency converter;
- Mini-type design, easy to clean, set and maintain;
- Can remove impurities cleanly.

Key Specification:

- | | |
|------------------------------|--------------------------|
| ● Control | Frequency converter |
| ● Length of machinable fiber | 22~42mm |
| ● Feed volume | About 0.25-1 kg per hour |
| ● Cylinder | Φ300mm, 1500rpm |
| ● Dimension | 1800×800×1300mm(L×W×H) |