

Oxygen Concentrator



OXYGEN GENERATOR in Kunshan Dynmed Medical Technology Co., Ltd (www.dynmed.cn)

1) Dynamic Group has three factories near Shanghai and one of them is specializing in Oxygen Concentrator. And therefore, Kunshan Dynmed Medical Technology Co., Ltd is established.

2) Series of Oxygen Generator

DO2 Series: **DO2-5AH, DO2-5AM, DO2-8AM, DO2-10AM**

3) Main characters of DYNMED Oxygen Concentrators:

- >> **R&D ability:** 2 teams for developing and assuring the quality of each machine.
- >> **Multiple oxygen flow choice:** 0~2 liters, 0~3 liters, 0~5 liters, 0~8 liters, 0~10 liters
- >> **Factory warehouse:** 32,000 m², 82 workers in 3 production lines.
- >> **AVG oxygen purity:** 95%
- >> **Low sound level:** less than 45dB (A)
- >> **Net weight:** only 15.9kg (for 5 liters), 19.2kg (for 8 liters), 22.6kg (for 10 liters)
- >> **Strong ABS shell:** Anti-aging, allows max 200kg stand upon the machine.
- >> **3 years warranty** for each model.

4) Visible Features

Design:

Compact, fashionable and economical design with brilliant performance

Inner Structure:

Sequential inner structure so as to easy maintenance

Multi-Functions:

- >> Time accumulated
- >> 4 alarm systems--- power off alarm, pressure circulation fault alarm and compressor fault alarm, low oxygen concentration alarm (optional function).
- >> Sterilization, dust-removing, oil free and smell removing.
- >> LCD screen display
- >> Super-low noise, energy saving, environmental protected
- >> Stable working time: **24 hours continuously**
- >> Provide **3 years guarantee** for quality defective.

Typical Features:

Dual oxygen output, offer average flow oxygen for two persons simultaneously. This can be customized for models-- **DO2-8AM and DO2-10AM.**

Oxygen humidifier:


Yes, all models are equipped with a humidifier bottle.

Optional functions:

- >> Oxygen sensor (for low oxygen concentration alarm)
- >> Nebulizer (allow inhaling oxygen treatment at the same time)

Easy operation:

You only need a nasal cannula or oxygen mask to connect the oxygen outlet and turn on the oxygen concentrator, then all is done! (Note: 1 pc nasal cannula is provided)





Kunshan Dynmed Medical Technology Co., Ltd

www.dynmed.cn

To create value for all our clients!

Note: This document is confidential, do not print or forward to a third party.

Dynmed Medical---Dynmed Medical---Dynmed Medical